



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

LED Surgical Lighting System

Powered by High Depth LEDs

Blitz
Celesta Series

Model:
Blitz Celesta-M



Specifications

- Light Intensity 200000 Lux
- Number of LEDs: 72
- Illumination Depth: 70mm
- Battery Backup
- Radiations: Zero IR and Zero UV radiations
- Colour Temperature: 4500 – 5000 K
- CRI: 95% - 98%
- Intensity control: 10% to 100%
- Temperature above Surgeon's Head: <1 Deg
- Height Adjustment: 1 meter
- LED Life: 60,000 Hours
- Sterilizable Handle provided



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad.

Ph: 9849466520, 7997192056

Email : lumensmedical@gmail.com www.lumensmedical.com

Distributor

ISO Certified, ISO 13485, CE, IEC 6060-1-1, IEC 6060-1-1-2, CDSCO CLASS A&B