



Germans call it  
Blitz  
We call it  
Pure Natural  
Surgical Light



Lumens Medical Pvt. Ltd.

*Life Enhancing Technology*

An ISO 9001, ISO - 13485, CE, WHO cGMP  
ISO -1400 Certified Company

## Blitz Series Hi Depth Surgical Illumination Systems

**POWERED BY HIGH DEPTH LED'S**

Surgical illumination systems with total shadow dilution from Lumens Medical Pvt Ltd. Lumens Medical Pvt Ltd, the most versatile LED Surgical lighting systems manufacturing company in the India brings in Blitz series LED Surgical illumination systems, with a Totally shadow free technology. The LED's life time of Blitz surgical lighting systems spans upto 20 years under standard usage.

*Perfect shadow dilution  
Natural Day light  
Laminar flow compatibility*

A Patent has Been filed by  
Lumens Medical Pvt Ltd  
For Blitz series  
illumination systems.

**Blitz series lights** Produces Hi depth & Low glare illumination with its unique combination of Yellow and white LED's resulting in a natural day light reproduction, Optimal depth of illumination achieved with no refocusing necessary, Provides Large illuminated area for most demanding surgical procedures.

- Long durability LED'S with 20 years of working life in standard conditions.
- Designed with smooth flowing contours for optimal laminar Air flow compatibility
- Removable sterilizable grip • Intensity control by Micro Steps. • Endo mode facility

## Blitz HD Single Dome ceiling light with 1,80,000 lux intensity



### Blitz HD - S

#### Specifications

No of LED'S Used	:	32
Lux Intensity	:	1,80,000 LUX
Intensity adjustment	:	soft touch 10-100%
Endo Mode	:	available
Colour temperature	:	4,500 – 5,200
CRI	:	>95ra
Display	:	Clear LCD Screen
Temperature increase over surgeons head	:	<2 Degrees
Index spot for surgery	:	200-250mm
Diameter of Light head	:	400mm each light head



**Lumens Medical Pvt. Ltd.**

Illumination Technology Division

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893

Email : lumensmedical@gmail.com www.lumensmedical.com

In technical collaboration with PHOTONIC-Austria

Area Distributor