



Lumens Medical Pvt. Ltd.™

Life Enhancing Technology

ISO Certified, ISO 13485,
CE, IEC, ISO 60601-2-41:2021

Blitz L6 Series

LED Surgical Lighting System

High Depth Super Bright LEDs



**Model:
Blitz L6 Twin Dome**

- Light Intensity 100000 Lux + 100000 Lux
- Battery Backup (Optional)
- Colour Temperature: 3500 – 5000 K adjustable
- Intensity control: 10% to 100%
- Height Adjustment: 1 meter
- Sterilizable Handle provided
- Illumination Depth: 70mm
- Radiations: Zero IR and Zero UV radiations
- CRI: 95%
- Temperature above Surgeon's Head: <math><1^{\circ}\text{C}</math>
- LED Life: 60,000 Hours
- Remote Control provided (Optional)

#2-24-113/3/B/NR, Laxmi Narayana Nagar, IDA, Uppal, Hyderabad-500039, TS, India.

Ph: 7997192056, 9849466520