



Lumens Medical Pvt. Ltd.

*Life Enhancing Technology*

An ISO 9001, ISO - 13485, CE, WHO cGMP  
ISO -1400 Certified Company

# *Blitz M Series*

## *Hi Depth Surgical Illumination Systems*

**POWERED BY HIGH DEPTH LED'S**

Surgical illumination systems with total shadow dilution from Lumens Medical Pvt Ltd. Lumens Medical Pvt Ltd, the most versatile LED Surgical lighting systems manufacturing company in the India brings in Blitz series LED Surgical illumination systems, with a Totally shadow free technology. The LED's life time of Blitz surgical lighting systems spans upto 20 years under standard usage.

**Germans call it Blitz - We call it Pure Natural Surgical Light**

*Perfect shadow dilution*  
*Natural Day light*  
*Laminar flow compatibility*

**PATENTED**

**with No. 308532**

by Controller General of Patents,  
Designs, and Trade Marks  
Ministry of Commerce & Industry

*Blitz series lights M6 Produces Hi depth & Low glare illumination with its unique combination of Yellow and white LED's resulting in a natural day light reproduction, Optimal depth of illumination achieved with no refocusing necessary, Provides Large illuminated area for most demanding surgical procedures.*

- Long durability LED'S with 20 years of working life in standard conditions.
- Designed with smooth flowing contours for optimal laminar Air flow compatibility
- Removable sterilizable grip • Brightness control built-in on light dome itself • Endo mode facility

## **Blitz Mobile light with 1,60,000 lux intensity**

**Model:**  
**Blitz ML6**



### **Specifications**

No of LED'S Used	: 54
Lux Intensity	: 1,60,000 LUX
Intensity adjustment	: soft touch 10-100%
Endo Mode	: available
Colour temperature	: 4,500 – 5,200
CRI	: >95 Ra
Display	: Clear LCD Screen
Temperature increase over surgeons head	: <2 Degrees
Index spot for surgery	: 200-250mm
Diameter of Light head	: 400mm each light head



**Lumens Medical Pvt. Ltd.**

Illumination Technology Division

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893

Email : lumensmedical@gmail.com www.lumensmedical.com

In technical collaboration with PHOTONIC-Austria

Area Distributor