

# Mobile Surgical Light

*Indigenously manufactured LM Series surgical OT lights with Laminar flow compatibility, compact size with pure white light and shadow free illumination.*

*The LM Series lights come with different lux intensities depending on surgeons requirement.*

**Model : M-120**



Quality Certificates  
ISO 9001 : 2015  
ISO 13485 : 2012  
ISO 1400 : 2015  
CE, WHO cGMP



## Specifications :

- *Lux intensity 1,20,000 Lux*
- *Colour Temperature 4,500- 5200K*
- *CRI >95 Ra*
- *Intensity adjustment 10% - 100% on LCD Screen Display*
- *Endo-mode*
- *Temperature increase over surgeons head not >2 Degrees.*
- *Index light spot 200mm-250mm*
- *Power 230V 50 Hz*
- *Dome diameter 400mm*



Photonic Optische Geräte GmbH & Co KG  
Seeböckgasse 59  
1160 Wien, Austria



**Lumens Medical Pvt. Ltd.**<sup>TM</sup>

Illumination Technology Division

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893

Email : lumensmedical@gmail.com www.lumensmedical.com

In technical collaboration with PHOTONIC-Austria