



Lumens Medical Pvt. Ltd.

# DOUBLE DOME SURGICAL OT LIGHT

Indigenously manufactured LM Series surgical OT lights with Laminar flow compatibility, compact size with pure white light and shadow free illumination. The LM Series lights come with different lux intensities depending on surgeons requirement.



Quality Certificates  
ISO 9001 : 2015  
ISO 13485 : 2012  
ISO 1400 : 2015  
CE, WHO cGMP

Model: LMDD 280

## Technical Specifications:

- Lux intensity 2,80,000 Lux
- Colour Temperature 4,500- 5200K
- CRI >95 Ra
- Intensity adjustment 10% - 100% on LCD Screen Display.
- Endo-mode
- Temperature increase over surgeons head not >2 Degrees.
- Index light spot 200mm-250mm
- Power 230V 50 Hz



Lumens Medical Pvt. Ltd.

Illumination Technology Division

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893  
Email : lumensmedical@gmail.com www.lumensmedical.com

Area Distributor