



**Lumens Medical Pvt. Ltd.**

*Life Enhancing Technology*

An ISO 9001, ISO - 13485, CE, WHO cGMP  
ISO -1400 Certified Company

**Model : LMSD -140**

***Designed for Optimum  
Laminar Air Flow Compatibility***



***Specifications :***

- Illumination Intensity: 1,40,000 Lux
- Height Adjustment: 600mm
- Action Radius: 1850mm
- Smooth Radial, Angular & Axial movements
- LED life: 50,000 Hrs
- Colour Rendering Index (CRI) : >95 Ra.
- Endo Mode on One Touch
- Intensity/ Brightness/ Contrast, Power switch, Endo mode control on Dome Handle
- Temperature raise over surgeons head <10 Deg
- Colour Temperature: 4,500-5,200 K
- Intensity Adjustment: 10% - 100% on LCD Screen Display
- Temperature Increase over Surgeon's Head Area: <2°C
- Index Light Spot: 200 mm - 250 mm
- Power: 230 VAC 50Hz



***Indigenously manufactured LM Series surgical OT lights with Laminar flow compatibility, compact size with pure white light and shadow free illumination. The LM Series lights come with different lux intensities depending on surgeons requirement.***



**Lumens Medical Pvt. Ltd.**

**Illumination Technology Division**

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893

Email : lumensmedical@gmail.com www.lumensmedical.com

**In technical collaboration with PHOTONIC-Austria**



ISO 9001:2008 NR.16459/0  
EN ISO 13485:2012 NR.00192/0

Photonic Optische Geräte GmbH & Co KG  
Seeböckgasse 59  
1160 Wien, Austria

