



Mobile Surgical Light

Indigenously manufactured LM Series surgical OT lights with Laminar flow compatibility, compact size with pure white light and shadow free illumination.

The LM Series lights come with different lux intensities depending on surgeons requirement.

Model : M-140



Quality Certificates
ISO 9001 : 2015
ISO 13485 : 2012
ISO 1400 : 2015
CE, WHO cGMP



Specifications :

- **Lux intensity 1,40,000 Lux**
- **Colour Temperature 4,500- 5200K**
- **CRI >95 Ra**
- **Intensity adjustment 10% - 100% on LCD Screen Display**
- **Endo-mode**
- **Temperature increase over surgeons head not >2 Degrees.**
- **Index light spot 200mm-250mm**
- **Power 230V 50 Hz**



Photonic Optische Geräte GmbH & Co KG
Seeböckgasse 59
1160 Wien, Austria



Lumens Medical Pvt. Ltd.

illumination Technology Division

2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,
Telangana, India. Tel.: + 91 7997192056, 6309995893

Email : lumensmedical@gmail.com www.lumensmedical.com

In technical collaboration with PHOTONIC-Austria