



# Examination Lamp

for Labour Room, ENT,  
Dermatology, Plastic Surgery,  
OP Rooms, Minor Surgery Room etc

Swan neck adjustable  
in any desired angle

Intensity  
Adjustment


Model: L-40

Height  
Adjustment

Handle  
to  
move

Instrument  
Placing  
Platform

Wheels  
with brakes



## Specifications :

- Lux intensity 40,000 Lux
- Intensity adjustment by Rotary Knob
- Height adjustment facility
- Mobiles on wheels
- Instrument placing platform
- Power 230V 50 Hz



Photonic Optische Geräte GmbH & Co KG  
Seeböckgasse 59  
1160 Wien, Austria

Lumens Medical Pvt. Ltd.  
Illumination Technology Division

# 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad-500 039,  
Telangana, India. Tel.: + 91 7997192056, 6309995893  
Email : lumensmedical@gmail.com www.lumensmedical.com  
In technical collaboration with PHOTONIC-Austria