

Focus Lamp for Examination

Model: L70



Technical Specifications:

- Source of Light: LED Light
- Minimum 40,000 Lux at a Working distance of 0.5m
- Height adjustment: ≤ 440 mm
- Radial and axial movement of the lamp
- Mobile floor mounting type with 150cm minimum height including light head.
- Made of high quality ABS / Equivalent material body
- Heat dissipation: Maintain normal temperature and heat be disbursed through a cooling mechanism.
- Mobility, Portability: Mobile
- Battery Operated : Optional
- Power requirement: Recharging Unit: Input Voltage – 220V-240V AC, 50Hz



Lumens Medical Pvt. Ltd.TM

Life Enhancing Technology

2-24-113/3/B/NR, Laxmi Narayana Nagar, Uppal, Hyderabad-500 039, Telangana.

Tel.: 7997192056, 9324001836

Email : lumensmedical@gmail.com www.lumensmedical.com