

Germans call it  
**Blitz**  
We call it  
Pure Natural  
Surgical Light



Lumens Medical Pvt. Ltd.™  
*Life Enhancing Technology*

## Blitz Series Hi Depth Surgical Illumination Systems

**POWERED BY HIGH DEPTH LED'S**

"**Blitz Glacea X2**" is a High depth Surgical illumination System with total shadow dilution. The low glare illumination with high depth enables Surgeon to do complicated & time consuming Surgeries without any Eye strain.



Model: **Blitz Glacea X**

### Key features:

- Focus adjustment available
- Colour Temperature 4500 - 5200 K
- Intensity Adjustment 10% - 100%
- CRI >95
- Display clear LCD Screen
- LED Life 60,000 Hours



### SPECIFICATIONS

- Number of LEDs 48
- Lux Intensity 2,00,000 LUX
- Intensity Adjustment 10%-100% Soft touch
- Endo Mode Available
- CRI >95
- Colour Temperature 4500 – 5200 K
- Display Clear LCD Screen
- Temperature raise over surgeons head <1 Degree
- Index spot for surgery 150 mm – 300mm adjustable
- LED Life 60,000 Hrs
- Illumination depth up to 70mm
- Battery backup (Optional)
- Height adjustment 1 meter
- Remote control (Optional)

*Our Pride*  
 Over 2000 surgeries  
 Performed under our  
 Surgical lights per day

*Installation base*  
 More than 1000 satisfied  
 surgeons using Lumens  
 Surgical lighting  
 systems all over India

**The Company:** Lumens Medical Private Limited, India's most versatile LED Surgical illumination systems manufacturer with more than 20 models of Operation theatre lighting systems. With its own in house R&D team.

**Quality Certifications**  
 ISO 9001  
 ISO 13485  
 IEC 60606-1-1  
 IEC 6060-1-1-2  
 CDSCO



**MADE IN INDIA  
 FOR THE WORLD**



**Lumens Medical Pvt. Ltd.**  
*Life Enhancing Technology*

An ISO Certified, ISO-13485, CE, IEC 60601-1, IEC 60601-1-2, CDSCO CLASS A&B  
 # 2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad.  
 Tel.: 9849466520, 7997192056  
 Email : lumensmedical@gmail.com www.lumensmedical.com