



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

LED Surgical Lighting System

Powered by High Depth LEDs

Blitz Series

Model:
Blitz Celesta-D



Laminar flow Compatibility
Colour Fusion Technology
VIBGYOR
Rainbow spectrum
for achieving Natural
Day light on OT Table


Colour Lights

Remote sensor

Intensity controllers

Removable, Sterilizable Handle

Yellow Lights



ISO Certified, ISO 13485, CE, IEC 6060-1-1, IEC 6060-1-2, CDSCO CLASS A&B

Surgical illumination systems with total shadow dilution from **Lumens Medical Pvt. Ltd.** The most versatile LED surgical lighting systems manufacturing company in India brings in Blitz Series LED Surgical illumination Systems with a totally shadow free technology. The LED's life time of Blitz Surgical Lighting Systems spans upto 20 years under standard usage.

*Perfect shadow dilution
Natural Day Light
Laminar flow compatibility*



Also available in
Yellow & White Light
mixing of colour

Specifications of Blitz Celesta-D

- Light Intensity 200000 Lux + 200000 Lux
- Number of LEDs: 72 + 72
- Illumination Depth: 70mm
- Battery Backup
- Radiations: Zero IR and Zero UV radiations
- Colour Temperature: 4500 – 5000 K
- CRI: 95% - 98%
- Intensity control: 10% to 100%
- Temperature above Surgeon's Head: <1 Deg
- Height Adjustment: 1 meter
- LED Life: 60,000 Hours
- Sterilizable Handle provided



Lumens Medical Pvt. Ltd.TM

Life Enhancing Technology

An ISO 9001, ISO - 13485, CE, CDSO CLASS A&B Certified Company

2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad.

Ph: 9849466520, 7997192056

Email : lumensmedical@gmail.com www.lumensmedical.com

Germans
call it Blitz
We call it
Pure Natural
Surgical
Light

Distributor