

Fully Electric C-Arm OT Table

For Orthopedic, Gynaec & General Surgeries



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology
An ISO 9001, ISO-13485, CE, ISO-1400 Certified Company

Model: **LM-OT1E**



Horse Shoe
Neuro Attachments



Ortho Attachments

Features:

- Manual Mechanism for functioning and controlling the table.
- Centric base made from high quality stainless steel.
- LEG and HEAD SECTION are detachable and interchangeable.
- Inbuilt Kidney Bridge.
- Base and cover made of attractive, easy to clean stainless steel.
- Sophisticated mechanism provides smooth, step less articulation of the table top for precise patient positioning.

Technical Specifications:

- All positions electrically operated
- Top Dimension: L 1900mm x W 520mm
- Height Adjustment: 762mm -1000 mm
- Trendelenburg / Reverse: 30° / 25°
- Lateral Tilt: 20° / 20°
- Kidney Elevator: 130mm
- Power Supply: 24V DC
- Back Rest (up / down): 80° / 25°
- Leg Rest (up / down): 10° / 90°
- Head Rest (up / down): 45° / 90°
- Battery Backup: 2-3 Hours (optional)
- Mains Voltage Capacity: AC 220V±10%50Hz
- Patient Weight Capacity: 200 Kgs



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

An ISO 9001, ISO - 13485, CE, WHO cGMP ISO -1400 Certified Company

2-24-113/3/B/NR, 1st Floor, Laxmi Narayana Nagar, Uppal, Hyderabad - 500039

Ph: 9849466520, 7997192056

